



Catalogue - Price List

Design & Production:

PolyDiagnost GmbH
Schäfflerstr. 13
D- 85276 Pfaffenhofen/Ilm
Germany
Phone +49(8441)8060-0
Fax +49(8441)8060-60

info@polydiagnost.com
www.polydiagnost.com

Worldwide Sales & Marketing:


AD-TECH-MED Ltd.
Orkana 23 / 43
PL-20504 Lublin
Poland
Phone +48(81)5322461
Fax +48(81)5322462

office@adtechmed.eu
www.adtechmed.eu








The PolyDiagnost rigid optic modular system combination table

The interdisciplinary rigid endoscope- one optic for every application!
Without any decontamination issues due to the disposable PolyShaft.





These are recommended combinations. For other combinations and fields of applications feel free to contact us.

Optics

Order No.	Article Description	Additional Description	Picture
PD-DS-1080	rigid optic modular 30,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application highest resolution examination on the highest level no decontamination issues 30,000 pixel / 70° view outer diameter: 1.3 mm standard light connection: ACMI	
PD-DS-1083	rigid optic modular 10,000 Pixel / 120°	the interdisciplinary rigid endoscope one optic for every application high resolution at low diameter, 120° wide angle no decontamination issues 10,000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
PD-DS-1084	rigid optic modular 6,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application very thin with best resolution for examination and therapy no decontamination issues 6,000 pixel / 70° view outer diameter: 0.55 mm standard light connection: ACMI	
PD-DS-1086	rigid optic modular 3,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application ultrathin, with finest resolution especially for diagnostic in thinnest channels no decontamination issues 3,000 pixel / 70° view outer diameter: 0.45 mm standard light connection: ACMI	
PD-PN-0080	rigid optic modular PCNL 10,000 Pixel / 120°	extra long rigid endoscope special design for PCNL high resolution at low diameter, 120° wide angle no decontamination issues 10,000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	












Optics

Order No.	Article Description	Additional Description	Picture
PD-PS-0095	flexible optic PolyScope 120° view / 300cm	PolyScope- the flexible optic system no decontamination issues due to disposable catheters 10,000 Pixel, 120° View, 300cm length	
PD-PS-0096	flexible optic PolyScope 120° view / 180cm	PolyScope- the flexible optic system no decontamination issues due to disposable catheters 10,000 Pixel, 120° View, 180cm length	

Accessories for Optics

Order No.	Article Description	Additional Description	Picture
PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm best for use with: 10,000 and 30,000 pixel optic	
PD-FS-4004	ocular fix-focus blue	modular ocular fix focus blue with pin for easy mounting at the 3-joint-arm best for use with: 6,000 and 3,000 pixel optic	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-0097	3-joint-arm PolyScope	for use with flexible optics consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting 2 pcs.	
PD-VC-0220	cable clip set	set of 10 pieces for fixation of cables to the 3-joint-arm	
PD-AC-1012	SteriTray large	with individual holders for the endoscope suitable for PCNL optic	
PD-AC-1013	SteriTray small	with individual holders for the endoscope suitable for rigid optics (except PD-PN-0080)	

Rigid Endoscopy


Order No.	Article Description	Additional Description	Picture
PD-DS-1315	handle middle	suitable for PolyShaft 50mm for TMJ, TearChannel, ECP curved	
PD-DS-1320	handle short	suitable for PolyShaft 90mm for Saliva Channel, Milk Duct, PNL	
PD-TD-0030	cleaning brush handle set of 5	cleaning brush with eye for all handles set of 5 pieces	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-CA-0610	Acha 6,000 Pixel	for use of rigid optic without PolyShaft suitable for 6.000 pixel optic diameter: 0.59mm LuerLock	
PD-CA-0611	Acha 10,000 Pixel	for use of rigid optic without PolyShaft for 10.000 pixel optic diameter: 0.95mm LuerLock	
PD-CA-0612	Acha 30,000 Pixel	for use of rigid optic without PolyShaft for 30.000 pixel optic diameter: 1.33mm LuerLock	
PD-DS-1000	PolyShaft straight 0.65 / 90 mm	Outer diameter: 0.65mm Working length: 90mm disposable / sterile packed Suitable for 3.000 pixel optic only on request	 similar to illustration
PD-DS-1005	PolyShaft straight 0.75 / 90 mm	Outer diameter: 0.75mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	 similar to illustration
PD-DS-1010	PolyShaft straight 0.90 / 90 mm	Outer diameter: 0.90mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	

Please notice: errors excepted; subject and prices to change without notice; similar illustrations












Rigid Endoscopy

Order No.	Article Description	Additional Description	Picture
PD-DS-1015	PolyShaft straight 1.15 / 90 mm	Outer diameter: 1.15mm Working length: 90mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1020	PolyShaft straight 1.60 / 90 mm	Outer diameter: 1.60mm Working length: 90mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1025	PolyShaft straight 2.00 / 90 mm	Outer diameter: 2.00mm Working length: 90mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-DS-1030	PolyShaft straight 0.75 / 50 mm	Outer diameter: 0.75mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-DS-1035	PolyShaft straight 0.90 / 50 mm	Outer diameter: 0.90mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1040	PolyShaft straight 1.15 / 50 mm	Outer diameter: 1.15mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1042	PolyShaft straight 16G 1.60 / 50 mm	TMJ / Athroscopy Outer diameter: 1.60mm (16G) Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-DS-1045	PolyShaft curved 0.75 / 50 mm	Outer diameter: 0.75mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-DS-1050	PolyShaft curved 0.90 / 50 mm	Outer diameter: 0.90mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-DS-1055	PolyShaft curved 1.15 / 50 mm	Outer diameter: 1.15mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration

Rigid Endoscopy


Order No.	Article Description	Additional Description	Picture
PD-DS-1065	PolyShaft straight 18G 1.30/ 50 mm sharp tip	TMJ / Arthrocentesis sharp tip Outer diameter: 1.30mm (18G) Working length: 50mm disposable / sterile packed packing unit: 100pcs	 similar to illustration
PD-DS-1185	Trocar sharp	for PolyShaft 16G autoclavable	
PD-DS-1195	Obturator blunt	for PolyShaft 16G autoclavable	
PD-AC-0051	Solex - nipple- expander 0.95 mm	essential for milk duct endoscopy complete with 50 obturators suitable for PolyShaft 0.9mm	
PD-AC-0050	Solex - nipple- expander 1.1 mm	essential for milk duct endoscopy complete with 50 obturators suitable for PolySchaft 1.15mm	
PD-AC-0052	Solex - nipple- expander 1.3 mm	essential for milk duct endoscopy complete with 50 Obturators	
PD-AC-0055	obturator	for Solex packing unit: 50pcs	
PD-AC-0060	breast marker rhombus	outer diameter: 0.35mm, working length: 400mm please order separately: insertion tool PD-AC-0062	
PD-AC-0061	breast marker hook	outer diameter: 0.35mm; working length: 400mm please order separately: insertion tool PD-AC-0062	
PD-AC-0062	insertion tool	necessary for: breast markers and cytology brushes PD-AC-0060, PD-AC-0061, PD-TI-0005	

Rigid Endoscopy








Order No.	Article Description	Additional Description	Picture
PD-PN-0015	All-Seeing Needle	consisting of: Needle, Working Shaft 4.8F, Mandrin disposable / sterile packed	
PD-PN-0005	8F Working Shaft	OD 2,66mm 8F ID 2,3mm 7F including Dilator autoclavable	
PD-PN-0020	10,5F Working Shaft	OD 3,46mm 10,5F ID 3,00mm 8,5F including Dilator autoclavable	
PD-PN-0010	Lunderquist guide wire	Size 0.035 x 80cm	
PD-PN-0030	2-port adaptor	2-port adaptor for PCNL - 10,5 F working shaft	
PD-PN-0007	complete adapter set PCNL	consisting of: a) T-piece b) 3-port with female Luer-Lock c) Luer-Lock cap d) Tuohy Borst e) perfusion tube	
PD-CA-0104	Y-adaptor with haemostasis valve	disposable / sterile packed	
PD-CA-0112	Tuohy Borst-Adaptor	disposable / sterile packed	
PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
PD-CA-0202	rubber sluice	disposable / sterile packed	
PD-CA-0203	Y-adaptor	disposable / sterile packed	

Please notice: errors excepted: subject and prices to change without notice: similar illustrations


Rigid Endoscopy

Order No.	Article Description	Additional Description	Picture
PD-CA-0205	T-piece	disposable / sterile packed	
PD-CA-0601	Y-adaptor with check valve	1 set consisting of: 1x PD-CA-0606 1x PD-CA-0203s disposable / sterile packed	
PD-CA-0604	micro PNL set	consisting of: Y- adaptor T- piece rubber sluice disposable / sterile packed	
PD-CA-0606	check valve	required for drill PD-DS-1620 disposable / sterile packed	
PD-CA-0607	3-Port female with rotating LuerLock male	disposable / sterile packed	
PD-CA-0609	LuerLock cap set	set 5 pieces disposable / sterile packed	
PD-DS-0500	lavage pipe 0.5	Lavagepipe incl. check valve outer diameter: 0.5mm	
PD-DS-0505	lavage pipe 0.8	Lavagepipe incl. check valve outer diameter: 0.8mm	
PD-OI-0001	dilator	according to Wilder	
PD-OI-0010	Bowman- Probe 0.7	diameter: 0.7 mm Fig. 00/0	

Rigid Endoscopy

Order No.	Article Description	Additional Description	Picture
PD-OI-0011	Bowman- Probe 0.8	diameter: 0.8 mm Fig. 0/1	
PD-OI-0012	Bowman- Probe 0.9	diameter: 0.9 mm Fig. 1/2	
PD-OI-0013	Bowman- Probe 1.1	diameter: 1.1 mm Fig. 2/3	
PD-PS-0124	Laser- and Instrument shifter	For lasers with outer diameter 500-1000 μ	
PS-PS-0125	Laser- and Instrument shifter	For lasers with outer diameter 200-500 μ	
PD-DS-5102	forceps with scissorhandle	forceps, flexible with scissorhandle outer diameter: 0.8mm, length: 28cm	
PD-DS-1620	spiral drill	for handle middle with PolyShaft 50mm diameter: 0.38 mm / length over all: 193.8mm only to use with check valve PD-CA-0606 (1 included) autoclavable	
PD-DS-1630	spiral drill	for handle middle with PolyShaft 90mm diameter: 0.38 mm / length over all: 193.8mm only to use with check valve PD-CA-0606 (1 included) autoclavable	

Flexible Endoscopy

Order No.	Article Description	Additional Description	Picture
PD-PS-0144p	Modular Catheter 8 F 185cm	for Optic 300 cm, working length: 185 cm integrated light, working channel 1.2mm disposable / sterile packed	
PD-PS-0153p	Modular Catheter 8 F 30cm	for Optic 150 cm, working length: 30 cm integrated light, working channel 1.2mm disposable / sterile packed	
PD-PS-0154p	Modular Catheter 8 F 42cm	for Optic 150 cm, working length: 42 cm integrated light, working channel 1.2mm disposable / sterile packed	
PD-PS-0155p	Modular Catheter 8 F 20cm	for Optic 150 cm, working length: 20 cm integrated light, working channel 1.2mm disposable / sterile packed	
PD-PS-0174p	Modular Catheter 8 F 70cm	for Optic 180 cm, working length: 70 cm integrated light, working channel 1.2mm disposable / sterile packed	
PD-LC-9551	adaptor separate light	with pin for easy mounting at the 3- joint- arm standard light connection: ACMI	
PD-PS-0120	optic- combishifter		
PD-PS-0124	laser- and instrument shifter	For lasers with outer diameter 500-1000µ	
PS-PS-0125	laser- and instrument shifter	For lasers with outer diameter 200-500µ	

Electronics

Order No.	Article Description	Additional Description	Picture
PD-TR-13953	PolyPort rack 230V	with medically approved power strip for 230 Volt	
PD-TR-1500	PolyPort rack 110V	analogue to PD-TR-13953 for 110 Volt	
PD-VR-0100	EndoGnost® RA 100	The new all- in- one compact system high resolution camera with high speed shutter, xenon light- source, brilliant 15" monitor and recording function	
PD-VC-0210	camera EndoGnost® CS 207	High performance CCD camera for all endoscopic applications advanced brightness control with high speed shutter and Moirée filter suitable for 230 Volt / PAL	
PD-VC-0205	camera EndoGnost® CS 207	analogue to PD-VC-0210 for 110 Volt / NTSC	
PD-VC-1002	TV- adaptor 35 mm with snap fixture	with automatic snap fixture	
PD-VC-1022	TV- adaptor 35 mm autoclavable	soakable and autoclavable for easy cleaning	
PD-VR-0401	Medicapture USB 200	Easy full- featured video recording solution saves videos and images from the endoscope to USB stick Video in/out: S- video or composite BNC standard 3.5mm mini jack for foot switch (not included) connection please order seperately: PD-VR-0402 (foot switch)	












Electronics

Order No.	Article Description	Additional Description	Picture
PD-VM-3015w	EndoGnost® TFT monitor 15"	medical approved 15" TFT monitor with glass front and solid metal foot weight: 4,8 kg input: S- Video, DVI, VGA, RCA jack (cinch) voltage: 100- 240V	
PD-LS-0220	EndoGnost® LS 200 xenon light source	Xenon light source for a bright, sun- like daylight (5800°K) with brightness fine adjustment xenon lamp with high light yield and life cycle suitable for 230 Volt with IR/UV/Temp. Filter Suitable for all applications incl Ophthalmology	
PD-LS-0190	EndoGnost® LS 200 xenon light source	analogue to PD-LS-0220 for 110 Volt	
PD-IP-0100	irrigation pump IP 200	100ml/min suitable for 230 Volt incl. 5 pcs. PD-IP-0120 and 5 pcs. PD-IP-0121	
PD-IP-0110	irrigation pump IP 200	analogue to PD-IP-0100 for 110 Volt	
PD-OA-5129	micro drill system complete system	TCM 3000 / 230V for: rigid optic modular consisting of drive assembly with foot switch, drilling shaft and holder for the 3- joint- arm	
PD-OA-5128	micro drill system complete system	analogue to PD-OA-5129 for 110 V	






Accessories for Electronics

Order No.	Article Description	Additional Description	Picture
PD-VC-1100	dust cap for TV- adaptor		
PD-VR-0402	footswitch for Medicapture		
PD-LC-2001	light cable	length: 180 cm fiber bundle 4.8 mm without Adaptor	
PD-LC-3001	light cable kit	consists of: Light Cable PD-LC-2001 Adaptor Storz Lightsource PD-LC-9001 Adaptor ACMI Endoscope PD-LC-9101	
PD-LC-9001	adaptor C-LS Storz	adaptor for light cable PD-LC-2001 to fit to a Storz light source connection	
PD-LC-9002	adaptor C-LS Fuji	adaptor for light cable PD-LC-2001 to fit to a Fuji light source connection	
PD-LC-9003	adaptor C-LS Pentax	adaptor for light cable PD-LC-2001 to fit to a Pentax light source connection	
PD-LC-9004	adaptor C-LS Olympus/ACMI	adaptor for light cable PD-LC-2001 to fit to a Olympus/ACMI light source connection	
PD-LC-9005	adaptor C-LS Wolf	adaptor for light cable PD-LC-2001 to fit to a Wolf light source connection	
PD-LC-9006	adaptor C-LS Stryker	adaptor for light cable PD-LC-2001 to fit to a Stryker light source connection	
PD-LC-9101	adaptor C-EC ACMI	adaptor for light cable PD-LC-2001 to fit to a ACMI endoscope connection	









Accessories for Electronics

Order No.	Article Description	Additional Description	Picture
PD-LC-9103	adaptor C-EC Wolf	adaptor for light cable PD-LC-2001 to fit to a Wolf endoscope connection	
PD-LC-9104	adaptor C-EC Storz	adaptor for light cable PD-LC-2001 to fit to a Storz endoscope connection	
PD-LC-9502	adaptor EC Wolf	adaptor for a PolyDiagnost optic to fit to a Wolf light cable	
PD-LC-9503	adaptor EC Storz	adaptor for a PolyDiagnost optic to fit to a Storz light cable	
PD-EC-1411	cable kit for camera - monitor	1 x BNC 2m 1 x S-video 2m 1 x Adaptor BNC/Cinch	
PD-EC-1450	S- video- cable	length: 2 m	
PD-EC-1452	extension cable S- video	length: 5 m	
PD-EC-1480	cable BNC-BNC 1m	length: 1 m	
PD-EC-1481	cable BNC-BNC 2m	length: 2 m	
PD-EC-1483	cable BNC-BNC 5m	length: 5 m	
PD-EC-1486	adaptor Cinch/BNC		

Accessories for Electronics

Order No.	Article Description	Additional Description	Picture
PD-EC-1487	adaptor BNC/Cinch		
PD-EC-1488	coupling BNC/BNC		
PD-EC-1489	T- piece BNC		
PD-IP-0120	tubes for irrigation pump	1 piece	
PD-IP-0121	spike for tube	1 piece	












Instruments

Order No.	Article Description	Additional Description	Picture
PD-IC-0005	VascoCath 5.5F/60cm	for efficient and fast working in the vascular system with two working channels (0.6mm) diameter: 5.5F working length: 60 cm	
		disposable / sterile packed	
PD-IC-0010	VascoCath 5.5F/110cm	for efficient and fast working in the vascular system with two working channels (0.6mm) diameter: 5.5F working length: 110 cm	
		disposable / sterile packed	
PD-IC-0020	VascoCath 8F/60cm	for efficient and fast working in the vascular system with two working channels (0.6 and 1.2mm) diameter: 8F working length: 60 cm	
		disposable / sterile packed	
PD-IC-0025	VascoCath 8F/110cm	for efficient and fast working in the vascular system with two working channels (0.6 and 1.2mm) diameter: 8F working length: 110 cm	
		disposable / sterile packed	
PD-IC-0050	VascoCath DS 5.5F	for efficient and fast working in the vascular system with double steering with two working channels (0.6mm) diameter: 5.5F working length: 110 cm	
		disposable / sterile packed	
PD-IC-0055	VascoCath DS 8F	for efficient and fast working in the vascular system with double steering with two working channels (0.6 and 1.2mm) diameter: 8F working length: 110 cm	
		disposable / sterile packed	
PD-IC-0100	guide wire VascoCath 0.028	1 piece diameter: 0.028" length: 260 cm sterile	
PD-IC-0110	guide wire VascoCath 0.018	1 piece diameter: 0.018" length: 260 cm sterile	











Instruments

Order No.	Article Description	Additional Description	Picture
PD-OI-0015	Bangerter Lacrimal Probe Cannula	0.8 x 60 mm	
PD-OI-0016	Juenemann Lacrimal Probe Cannula	1.0mm with handle	
PD-TD-0100	guide wire 0.035"	length: 450 cm	
PD-TD-0110	guide wire 0.018"	length: 450 cm	
PD-TD-0130	guide wire 0.022"	length: 450 cm	
PD-TD-0150	guide wire 0.014"	length: 260 cm	
PD-TD-0160	guide wire 0.016"	length: 260 cm	
PD-TD-0601	laser shifter	for Storz / Walz NOT for PolyScope Diameter 0.9mm	
PD-TD-0602	gasket tube for laser shifter 0.45	gasket tube for laser shifter PD-TD-0601 for fibers up to 0.45 mm	
PD-TD-0603	gasket tube for laser shifter 0.85	gasket tube for laser shifter PD-TD-0601 for fibers up to 0.85 mm	
PD-TI-0005	cytology brush	for milk duct and saliva channel 1 set consisting of 5 brushes and 5 rotators STERILE please order separately: insertion tool PD-AC-0062	

Instruments

Order No.	Article Description	Additional Description	Picture
PD-TI-0082	modular cytology brush complete	for PolyScope (epiduroscopy) consisting of: PD-TI-2106, PD-TI-0084, PD-TI-0085	
PD-TI-0084	overcover for modular cytology brush	length: 80 cm autoclavable	
PD-TI-0085	modular cytology brush	length: 100 cm disposable / sterile packed	
PD-TI-2002	cytology brush	diameter: 1.0 mm length: 120 cm compact version (with handle) for PolyScope 20, 30, 42, 70 cm	
PD-TI-0102	cytology brush for ERCP	diameter: 1.0 mm length: 220 cm compact version (with handle)	
PD-TI-0103	basket for ERCP	5 wire basket with diameter 10mm diameter: 1.0 mm length: 220 cm compact version (with handle)	
PD-TI-0104	modular basket for ERCP	4 wire basket with diameter 15mm diameter: 0.8 mm length: 250 cm please order separately: PD-TI-2106, PD-TI-0105, PD-TI-0106	
PD-TI-0105	overtube for modular basket ERCP	with Luer Lock disposable / sterile packed	
PD-TI-0106	insertion tool for modular basket	to insert basket into overtube suitable for: PD-TI-0104, PD-TI-2006	
PD-TI-0107	rotator 1mm for basket	suitable for: PD-TI-0104, PD-TI-2006 disposable	
PD-TI-2003	basket for Uretero, PTC	diameter: 1.0 mm length: 120 cm compact version (with handle) for Polyscope 20, 30, 42, 70 cm reusable, nonsterile	

Instruments

Order No.	Article Description	Additional Description	Picture
PD-TI-2006	modular basket for Uretero, PTC	diameter: 0.8 mm length: 125 cm please order seperately: PD-TI-2106, PD-TI-2008, PD-TI-0106	
PD-TI-2008	overtube for modular basket Uretero	with Luer Lock disposable / sterile packed	
PD-TI-2012	tipless basket for Urology	compact version incl. overtube and handle diameter: 0.59 mm length: 120 cm disposable / sterile packed	
PD-TI-2104	basket 3 wire	3 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	
PD-TI-2105	basket 4 wire	4 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	
PD-TI-2106	handle for modular basket with chuck	suitable for: PD-TI-2104, PD-TI-2105, PD-TI-0104, PD-TI-2006 autoclavable	
PD-TI-0100	biopsy forceps for ERCP 0.8	diameter: 0.8 mm length: 220 cm reusable, nonsterile	
PD-TI-0101	biopsy forceps for ERCP 1.0	diameter: 1.0 mm length: 220 cm reusable, nonsterile	
PD-TI-2001	biopsy forceps for PolyScope	diameter: 1.0 mm length: 120 cm compact version (with handle) for PolyScope 20, 30, 42, 70 cm reusable, nonsterile	
PD-TI-2004	stone grasping forceps for Uretero, PTC, Broncho	diameter: 1.0 mm length: 120 cm compact version (with handle) reusable, nonsterile	
PD-TI-4100	biopsy forceps 60cm	diameter: 0.8 mm length: 60 cm oval spoon suitable for rigid endoscopy reusable, nonsterile	

PD-TI-4101

biopsy forceps 30cm

diameter: 0.8 mm

length: 30 cm



oval spoon

suitable for rigid endoscopy











reusable, nonsterile



Instruments

Order No.	Article Description	Additional Description	Picture
PD-TI-4104	biopsy forceps flexible	diameter: 0.8 mm length: 60 cm round spoon reusable, nonsterile	
PD-TI-5102	grasping forceps 30cm	diameter:0.8 mm length: 30cm suitable for rigid endoscopy reusable, nonsterile	








Spare Parts

Order No.	Article Description	Additional Description	Picture
PD-AC-0101	table clamp		
PD-AC-0102	vario arm		
PD-AC-0103	camera support	camera support for 3-joint-arm with pins to attach mounting	
PD-AC-0114	mounting	to attach modular oculars and adaptor separate light	
PD-LS-9901	bulb for EndoGnost® LS 200	SPARE PART Type XBO R 100 W/45 OFR	
PD-IP-0122	foot switch for IP 200	SPARE PART	
PD-IP-0123	bottle for irrigation pump	SPARE PART	
PD-OA-5141	drill shifter		
PD-OA-5142	driveshaft		
PD-OA-5143	gear transmission piece		
PD-OA-5144	motor mount for 3-Joint-Arm		

Spare Parts for expired products

Order No.	Article Description	Additional Description	Picture
PD-OA-5214	spiral drill for NeurOptik T 6 (III.Gen.)	diameter 0,38/ L=317	
PD-OA-5216	spiral drill for SalivaScope (III.Gen.)	diameter 0,38 /L=227	
PD-OA-5220	spiral drill f. Vitr.T6short (III.G.)	diameter 0,38 / L=122	
PD-ZS-0615	Drill Salivascope Modular	Diameter 0.38 mm Length: 162 mm	
PD-ZS-0620	Drill VitrOptic Modular	Diameter 0.38 mm Length: 142 mm	
PD-ZS-0621	Drill VitrOptic Modular	Diameter 0.38 mm Length: 196 mm	
PD-TI-0110	Viewing Biopsy Forceps complete	for: Salivascope flex Laduscope T flex consisting of: PD-TI-0111, PD-TI-0113	
PD-ZS-0422	biopsy forceps integrated in cannula	Outer diameter: 1.1 mm Working length: 85 mm incl. handle suitable for rigid 6,000 pixel optic	
PD-ZS-0522	grasping forceps integrated in cannula	Outer diameter: 1.1 mm Working length: 85 mm incl. handle suitable for rigid 6,000 pixel optic	
PD-VR-0910	Spare Lamp for EndoGnost® LSCR 200		

Optics

Order No.	Article Description	Additional Description	Picture
PD-DS-1080	rigid optic modular 30,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application highest resolution examination on the highest level no decontamination issues 30,000 pixel / 70° view outer diameter: 1.3 mm standard light connection: ACMI	
PD-DS-1083	rigid optic modular 10,000 Pixel / 120°	the interdisciplinary rigid endoscope one optic for every application high resolution at low diameter, 120° wide angle no decontamination issues 10,000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
PD-DS-1084	rigid optic modular 6,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application very thin with best resolution for examination and therapy no decontamination issues 6,000 pixel / 70° view outer diameter: 0.55 mm standard light connection: ACMI	
PD-DS-1086	rigid optic modular 3,000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application ultrathin, with finest resolution especially for diagnostic in thinnest channels no decontamination issues 3,000 pixel / 70° view outer diameter: 0.45 mm standard light connection: ACMI	
PD-PN-0080	rigid optic modular PCNL 10,000 Pixel / 120°	extra long rigid endoscope special design for PCNL high resolution at low diameter, 120° wide angle no decontamination issues 10,000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
PD-PS-0095	flexible optic PolyScope 120° view / 300cm	PolyScope- the flexible optic system no decontamination issues due to disposable catheters 10,000 Pixel, 120° View, 300cm length	
PD-PS-0096	flexible optic PolyScope 120° view / 180cm	PolyScope- the flexible optic system no decontamination issues due to disposable catheters 10,000 Pixel, 120° View, 180cm length	

Accessories for Optics

Order No.	Article Description	Additional Description	Picture
PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm best for use with: 10,000 and 30,000 pixel optic	
PD-FS-4004	ocular fix-focus blue	modular ocular fix focus blue with pin for easy mounting at the 3-joint-arm best for use with: 6,000 and 3,000 pixel optic	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-0097	3-joint-arm PolyScope	for use with flexible optics consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting 2 pcs.	
PD-VC-0220	cable clip set	set of 10 pieces for fixation of cables to the 3-joint-arm	
PD-AC-1012	SteriTray large	with individual holders for the endoscope suitable for PCNL optic	
PD-AC-1013	SteriTray small	with individual holders for the endoscope suitable for rigid optics (except PD-PN-0080)	

Electronics

Order No.	Article Description	Additional Description	Picture
PD-TR-13953	PolyPort rack 230V	with medically approved power strip for 230 Volt	
PD-TR-1500	PolyPort rack 110V	analogue to PD-TR-13953 for 110 Volt	
PD-VR-0100	EndoGnost® RA 100	The new all- in- one compact system high resolution camera with high speed shutter, xenon light- source, brilliant 15" monitor and recording function	
PD-VC-0210	camera EndoGnost® CS 207	High performance CCD camera for all endoscopic applications advanced brightness control with high speed shutter and Moirée filter suitable for 230 Volt / PAL	
PD-VC-0205	camera EndoGnost® CS 207	analogue to PD-VC-0210 for 110 Volt / NTSC	
PD-VC-1002	TV- adaptor 35 mm with snap fixture	with automatic snap fixture	
PD-VC-1022	TV- adaptor 35 mm autoclavable	soakable and autoclavable for easy cleaning	
PD-VR-0401	Medicapture USB 200	Easy full- featured video recording solution saves videos and images from the endoscope to USB stick Video in/out: S- video or composite BNC standard 3.5mm mini jack for foot switch (not included) connection please order separately: PD-VR-0402 (foot switch)	

PD-VM-3015w EndoGnost® TFT monitor 15" medical approved 15" TFT monitor with glass front and solid metal foot
weight: 4,8 kg
input: S- Video, DVI, VGA, RCA jack (cinch)
voltage: 100- 240V



PD-LS-0220 EndoGnost® LS 200 xenon light source Xenon light source for a bright, sun- like daylight (5800°K) with brightness fine adjustment
xenon lamp with high light yield and life cycle suitable for 230 Volt
with IR/UV/Temp. Filter



Suitable for all applications incl Ophthalmology

PD-LS-0190 EndoGnost® LS 200 xenon light source analogue to PD-LS-0220 for 110 Volt

PD-IP-0100 irrigation pump IP 200 100ml/min suitable for 230 Volt
incl. 5 pcs. PD-IP-0120 and 5 pcs. PD-IP-0121





PD-IP-0110 irrigation pump IP 200 analogue to PD-IP-0100 for 110 Volt

PD-OA-5129 micro drill system complete system TCM 3000 / 230V for: rigid optic modular
consisting of drive assembly with foot switch, drilling shaft and holder for the 3- joint- arm














PD-OA-5128 micro drill system complete system analogue to PD-OA-5129 for 110 V

Accessories for Electronics

Order No.	Article Description	Additional Description	Picture
PD-VC-1100	dust cap for TV- adaptor		
PD-VR-0402	footswitch for Mediacapture		

PD-LC-2001	light cable	length: 180 cm fiber bundle 4.8 mm without Adaptor	
PD-LC-3001	light cable kit	consists of: Light Cable PD-LC-2001 Adaptor Storz Lightsource PD-LC-9001 Adaptor ACMI Endoscope PD-LC-9101	
PD-LC-9001	adaptor C-LS Storz	adaptor for light cable PD-LC-2001 to fit to a Storz light source connection	
PD-LC-9002	adaptor C-LS Fuji	adaptor for light cable PD-LC-2001 to fit to a Fuji light source connection	
PD-LC-9003	adaptor C-LS Pentax	adaptor for light cable PD-LC-2001 to fit to a Pentax light source connection	
PD-LC-9004	adaptor C-LS Olympus/ACMI	adaptor for light cable PD-LC-2001 to fit to a Olympus/ACMI light source connection	
PD-LC-9005	adaptor C-LS Wolf	adaptor for light cable PD-LC-2001 to fit to a Wolf light source connection	
PD-LC-9006	adaptor C-LS Stryker	adaptor for light cable PD-LC-2001 to fit to a Stryker light source connection	
PD-LC-9101	adaptor C-EC ACMI	adaptor for light cable PD-LC-2001 to fit to a ACMI endoscope connection	
PD-LC-9103	adaptor C-EC Wolf	adaptor for light cable PD-LC-2001 to fit to a Wolf endoscope connection	
PD-LC-9104	adaptor C-EC Storz	adaptor for light cable PD-LC-2001 to fit to a Storz endoscope connection	

PD-LC-9502	adaptor EC Wolf	adaptor for a PolyDiagnost optic to fit to a Wolf light cable	
PD-LC-9503	adaptor EC Storz	adaptor for a PolyDiagnost optic to fit to a Storz light cable	
PD-EC-1411	cable kit for camera - monitor	1 x BNC 2m 1 x S-video 2m 1 x Adaptor BNC/Cinch	
PD-EC-1450	S- video- cable	length: 2 m	
PD-EC-1452	extension cable S- video	length: 5 m	
PD-EC-1480	cable BNC-BNC 1m	length: 1 m	
PD-EC-1481	cable BNC-BNC 2m	length: 2 m	
PD-EC-1483	cable BNC-BNC 5m	length: 5 m	
PD-EC-1486	adaptor Cinch/BNC		
PD-EC-1487	adaptor BNC/Cinch		
PD-EC-1488	coupling BNC/BNC		

PD-EC-1489

T- piece BNC



PD-IP-0120

tubes for irrigation pump 1 piece



PD-IP-0121

spike for tube






1 piece




TMJ











Order No.	Article Description	Additional Description	Picture
PD-DS-1083	rigid optic modular 10.000 Pixel / 120°	the interdisciplinary rigid endoscope one optic for every application high resolution at low diameter, 120° wide angle no decontamination issues 10.000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-1013	SteriTray small	with individual holders for the endoscope suitable for rigid optics (except PD-PN-0080)	
PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm best for use with: 10.000 and 30.000 pixel optic	
PD-DS-1315	handle middle	suitable for PolyShaft 50mm for TMJ, TearChannel, ECP curved	
PD-TD-0030	cleaning brush handle set of 5	cleaning brush with eye for all handles set of 5 pieces	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-DS-1042	PolyShaft straight 16G 1.60 / 50 mm	TMJ / Athroscopy Outer diameter: 1.60mm (16G) Working length: 50mm disposable / sterile packed packing unit: 100pcs	 similar to illustration

Catalogue

PD-DS-1065	PolyShaft straight 18G 1.30/ 50 mm sharp tip	TMJ / Arthrocentesis sharp tip Outer diameter: 1.30mm (18G) Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-DS-1185	Trocar sharp	for PolyShaft 16G autoclavable	
PD-DS-1195	Obturator blunt	for PolyShaft 16G autoclavable	
PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
PD-CA-0609	LuerLock cap set	set 5 pieces disposable / sterile packed	












Tear Channel

Order No.	Article Description	Additional Description	Picture
PD-DS-1084	rigid optic modular 6.000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application very thin with best resolution for examination and therapy no decontamination issues 6.000 pixel / 70° view outer diameter: 0.55 mm standard light connection: ACMI	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-1013	SteriTray small	with individual holders for the endoscope suitable for rigid optics (except PD-PN-0080)	
PD-FS-4004	ocular fix-focus blue	modular ocular fix focus blue with pin for easy mounting at the 3-joint-arm best for use with: 6.000 and 3.000 pixel optic	
PD-DS-1315	handle middle	suitable for PolyShaft 50mm for TMJ, TearChannel, ECP curved	
PD-TD-0030	cleaning brush handle set of 5	cleaning brush with eye for all handles set of 5 pieces	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-DS-1030	PolyShaft straight 0.75 / 50 mm	Outer diameter: 0.75mm Working length: 50mm disposable / sterile packed packing unit: 100pcs	 similar to illustration

PD-DS-1035	PolyShaft straight 0.90 / 50 mm	Outer diameter: 0.90mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1040	PolyShaft straight 1.15 / 50 mm	Outer diameter: 1.15mm Working length: 50mm disposable / sterile packed	
		packing unit: 100pcs	
PD-TI-2104	basket 3 wire	3 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	 similar to illustration
PD-TI-2105	basket 4 wire	4 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	 similar to illustration
PD-TI-2106	handle for modular basket with chuck	suitable for: PD-TI-2104, PD-TI-2105, PD-TI-0104, PD-TI-2006 autoclavable	
PD-OA-5129	micro drill system complete system	TCM 3000 / 230V for: rigid optic modular consisting of drive assembly with foot switch, drilling shaft and holder for the 3- joint- arm	
PD-DS-1620	spiral drill	for handle middle with PolyShaft 50mm diameter: 0.38 mm / length over all: 193.8mm only to use with check valve PD-CA-0606 (1 included) autoclavable	
PD-CA-0606	check valve	required for drill PD-DS-1620 disposable / sterile packed	
PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
PD-CA-0609	LuerLock cap set	set 5 pieces disposable / sterile packed	

Saliva Channel

Order No.	Article Description	Additional Description	Picture
PD-DS-1084	rigid optic modular 6.000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application very thin with best resolution for examination and therapy no decontamination issues 6.000 pixel / 70° view outer diameter: 0.55 mm standard light connection: ACMI	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-1013	SteriTray small	with individual holders for the endoscope suitable for rigid optics (except PD-PN-0080)	
PD-FS-4004	ocular fix-focus blue	modular ocular fix focus blue with pin for easy mounting at the 3-joint-arm best for use with: 6.000 and 3.000 pixel optic	
PD-DS-1320	handle short	suitable for PolyShaft 90mm for Saliva Channel, Milk Duct, PNL	
PD-TD-0030	cleaning brush handle set of 5	cleaning brush with eye for all handles set of 5 pieces	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-DS-1005	PolyShaft straight 0.75 / 90 mm	Outer diameter: 0.75mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	similar to illustration
PD-DS-1010	PolyShaft straight 0.90 / 90 mm	Outer diameter: 0.90mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	

PD-DS-1015	PolyShaft straight 1.15 / 90 mm	Outer diameter: 1.15mm Working length: 90mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1020	PolyShaft straight 1.60 / 90 mm	Outer diameter: 1.60mm Working length: 90mm disposable / sterile packed	
		packing unit: 100pcs	
PD-DS-1025	PolyShaft straight 2.00 / 90 mm	Outer diameter: 2.00mm Working length: 90mm disposable / sterile packed	
		packing unit: 100pcs	similar to illustration
PD-AC-0050	Solex - nipple- expander 1.1 mm	essential for milk duct endoscopy complete with 50 obturators suitable for PolySchaft 1.15mm	
PD-AC-0051	Solex - nipple- expander 0.95 mm	essential for milk duct endoscopy complete with 50 obturators suitable for PolyShaft 0.9mm	
PD-OI-0001	dilator	according to Wilder	
PD-CA-0202	rubber sluice	disposable / sterile packed; Minimum order: 100 pcs.	
PD-TI-0005	cytology brush	for milk duct and saliva channel 1 set consisting of 5 brushes and 5 rotators STERILE please order separately: insertion tool PD-AC-0062	
PD-AC-0062	insertion tool	necessary for: breast markers and cytology brushes PD-AC-0060, PD-AC-0061, PD-TI-0005	
PD-TI-2104	basket 3 wire	3 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	 similar to illustration
PD-TI-2105	basket 4 wire	4 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	 similar to illustration
PD-TI-2106	handle for modular basket with chuck	suitable for: PD-TI-2104, PD-TI-2105, PD-TI-0104, PD-TI-2006 autoclavable	

PD-OI-0011	Bowman- Probe 0.8	diameter: 0.8 mm Fig. 0/1	
PD-OI-0012	Bowman- Probe 0.9	diameter: 0.9 mm Fig. 1/2	
PD-OI-0013	Bowman- Probe 1.1	diameter: 1.1 mm Fig. 2/3	
PD-TI-0107	rotator 1mm for basket	suitable for: PD-TI-0104, PD-TI-2006 disposable	
PD-TI-5102	grasping forceps 30cm	diameter:0.8 mm length: 30cm suitable for rigid endoscopy reusable, nonsterile	
PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
PD-CA-0609	LuerLock cap set	set 5 pieces disposable / sterile packed	
PD-OA-5129	micro drill system complete system	TCM 3000 / 230V for: rigid optic modular consisting of drive assembly with foot switch, drilling shaft and holder for the 3- joint- arm	
PD-OA-5216	spiral drill for SalivaScope (III.Gen.)	diameter 0,38 /L=227	
Recommended for better diagnostic			
PD-DS-1083	rigid optic modular 10.000 Pixel / 120°	the interdisciplinary rigid endoscope one optic for every application high resolution at low diameter, 120° wide angle no decontamination issues 10.000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
PD-FS-4005	Zoom Ocular (red)	modular ocular zoom ocular with pin for easy mounting at the 3-joint-arm best for use with: 6,000 and 3,000 pixel optic	

PD-FS-4001 ocular fix-focus black

modular ocular fix focus black with pin for easy mounting at the 3-joint-arm
best for use with:
10.000 and 30.000 pixel optic



Recommended for tiny channels

PD-DS-1086 rigid optic modular
3.000 Pixel / 70°

**the interdisciplinary rigid endoscope
one optic for every application**
ultrathin, with finest resolution
especially for diagnostic in thinnest channels
no decontamination issues
3.000 pixel / 70° view
outer diameter: 0.45 mm
standard light connection: ACMI



PD-DS-1000 PolyShaft straight
0.65 / 90 mm

Outer diameter: 0.65mm
Working length: 90mm
disposable / sterile packed


on request


similar to illustration


Milk Duct


Order No.	Article Description	Additional Description	Picture
PD-DS-1084	rigid optic modular 6.000 Pixel / 70°	the interdisciplinary rigid endoscope one optic for every application very thin with best resolution for examination and therapy no decontamination issues 6.000 pixel / 70° view outer diameter: 0.55 mm	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-1013	SteriTray small	with individual holders for the endoscope suitable for rigid optics (except PD-PN-0080)	
PD-FS-4004	ocular fix-focus blue	modular ocular fix focus blue with pin for easy mounting at the 3-joint-arm best for use with: 6.000 and 3.000 pixel optic	
PD-DS-1320	handle short	suitable for PolyShaft 90mm for Saliva Channel, Milk Duct, PNL	
PD-TD-0030	cleaning brush handle set of 5	cleaning brush with eye for all handles set of 5 pieces	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-DS-1010	PolyShaft straight 0.90 / 90 mm	Outer diameter: 0.90mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	
PD-DS-1015	PolyShaft straight 1.15 / 90 mm	Outer diameter: 1.15mm Working length: 90mm disposable / sterile packed packing unit: 100pcs	
PD-AC-0050	Solex - nipple- expander 1.1 mm	essential for milk duct endoscopy complete with 50 obturators suitable for PolySchaft 1.15mm	


Please notice: errors excepted; subject and prices to change without notice; similar illustrations


PD-AC-0051	Solex - nipple-expander 0.95 mm	essential for milk duct endoscopy complete with 50 obturators suitable for PolyShaft 0.9mm	
------------	------------------------------------	--	---


PD-AC-0060	breast marker rhombus	outer diameter: 0.35mm, working length: 400mm please order separately: insertion tool PD-AC-0062	
------------	--------------------------	--	---


PD-AC-0061	breast marker hook	outer diameter: 0.35mm; working length: 400mm please order separately: insertion tool PD-AC-0062	
------------	--------------------	--	---


PD-TI-0005	cytology brush	for milk duct and saliva channel 1 set consisting of 5 brushes and 5 rotators STERILE please order separately: insertion tool PD-AC-0062	
------------	----------------	---	---


PD-AC-0062	insertion tool	necessary for: breast markers and cytology brushes PD-AC-0060, PD-AC-0061, PD-TI-0005	
------------	----------------	---	---

PD-TI-2104	basket 3 wire	3 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	 similar to illustration
------------	---------------	---	---


PD-TI-2105	basket 4 wire	4 wire basket without handle suitable for rigid endoscopy autoclavable please order separately: handle for basket PD-TI-2106	 similar to illustration
------------	---------------	---	--

PD-TI-2106	handle for modular basket with chuck	suitable for: PD-TI-2104, PD-TI-2105, PD-TI-0104, PD-TI-2006 autoclavable	
------------	---	---	---

PD-CA-0201	perfusion tube	100cm with Luer Lock disposable / sterile packed	
------------	----------------	--	---

PD-CA-0609	LuerLock cap set	set 5 pieces disposable / sterile packed	
------------	------------------	--	---

Recommended for better diagnostic

PD-DS-1083	rigid optic modular 10.000 Pixel / 120°	the interdisciplinary rigid endoscope one optic for every application high resolution at low diameter, 120° wide angle no decontamination issues 10.000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
------------	--	--	---

PD-FS-4001 ocular fix-focus black modular ocular fix focus black with pin for easy mounting at the 3-joint-arm
best for use with:
10.000 and 30.000 pixel optic



Recommended for tiny breasts

PD-DS-1086 rigid optic modular 3.000 Pixel / 70° **the interdisciplinary rigid endoscope one optic for every application**
ultrathin, with finest resolution especially for diagnostic in thinnest channels
no decontamination issues
3.000 pixel / 70° view
outer diameter: 0.45 mm
standard light connection: ACMI



PD-DS-1000 PolyShaft straight 0.65 / 90 mm Outer diameter: 0.65mm
Working length: 90mm
disposable / sterile packed



on request

similar to illustration

PD-DS-1005 PolyShaft straight 0.75 / 90 mm Outer diameter: 0.75mm
Working length: 90mm
disposable / sterile packed



packing unit: 100pcs

similar to illustration

Recommended additional accessories

PD-LS-0220 EndoGnost® LS 200 xenon light source Xenon light source 100 W for a bright, sun-like daylight (5800°K) with brightness fine adjustment
xenon lamp with high light yield and life cycle suitable for 230 Volt
with IR/UV/Temp. Filter



Suitable for all applications incl Ophthalmology

PD-LC-3001 light cable kit consists of:
Light Cable PD-LC-2001
Adaptor Storz Lightsource PD-LC-9001
Adaptor ACMI Endoscope PD-LC-9101



Micro Perc PCNL

Order No.	Article Description	Additional Description	Picture
PD-PN-0080	rigid optic modular PCNL 10.000 Pixel / 120°	extra long rigid endoscope special design for PCNL high resolution at low diameter, 120° wide angle no decontamination issues 10.000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-AC-1012	SteriTray large	with individual holders for the endoscope suitable for PCNL optic	
PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm best for use with: 10.000 and 30.000 pixel optic	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-PN-0015	All-Seeing Needle	consisting of: Needle, Working Shaft 4.8F, Mandrin disposable / sterile packed	
PD-PN-0005	8F Working Shaft	OD 2,66mm 8F ID 2,3mm 7F including Dilator autoclavable	
PD-PN-0020	10,5F Working Shaft	OD 3,46mm 10,5F ID 3,00mm 8,5F including Dilator autoclavable	

PD-PN-0010 Lunderquist guide wire Size 0.035 x 80cm



PD-PN-0007 complete adapter set PCNL consisting of:
 a) T-piece
 b) 3-port with female Luer-Lock
 c) Luer-Lock cap
 d) Tuohy Borst
 e) perfusion tube



PD-IP-0100 Irrigation Pump IP 200 100ml/min
 220 Volt
 incl. 5 pcs. PD-IP-0120 and 5 pcs. PD-IP-0121



PD-IP-0120 Tubes for Irrigation Pump 1 piece



PD-IP-0121 Spike for Tube 1 piece

