







Optics

Order No.	Article Description	Additional Description	Picture
PD-PN-0080	rigid optic modular PCNL 10,000 Pixel / 120°	extra long rigid endoscope special design for PCNL high resolution at low diameter, 120° wide angle no decontamination issues 10,000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI AUTOCLAVABLE	
PD-AC-1012	SteriTray large	with individual holders for the endoscope suitable for PCNL optic	

Accessories for Optics

Order No.	Article Description	Additional Description	Picture
PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm best for use with: 10,000 and 30,000 pixel optic	
PD-AC-0096	3-joint-arm basic version	for use with rigid endoscopes consisting of: PD-AC-0101 table clamp PD-AC-0102 vario-arm PD-AC-0103 camera support PD-AC-0114 mounting	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	

PD-PN-1008 MicroPerc Set

MicroPerc set - sterile content: 3-port with female Luer-Lock, Luer-Lock cap, Tuohy Borst adaptor, perfusion tube, 3-way stopcock and 8F working shaft, All-seeing-needle, working shaft 4.8 F
working shaft 8F for the PCNL system suitable for the lithoclast® probe inner/outer diameter 7F/8F including Dilator




PD-PN-0020 10,5F Working Shaft OD 3,46mm 10,5F
ID 3,00mm 8,5F
including Dilator



autoclavable

PD-CA-0113 T-piece

disposable / sterile packed



PD-PN-0030 2-port adaptor 2-port adaptor for PCNL - 10,5 F working shaft
autoclavable



PD-LC-9503 adaptor EC Storz adaptor for a PolyDiagnost optic to fit to a Storz light cable

