













Micro Perc PCNL

Order No.	Article Description	Additional Description	Picture
PD-PN-0080	rigid optic modular PCNL 10.000 Pixel / 120°	extra long rigid endoscope special design for PCNL high resolution at low diameter, 120° wide angle no decontamination issues 10.000 pixel / 120° view outer diameter: 0.9 mm standard light connection: ACMI AUTOCLAVABLE	
PD-AC-1012	SteriTray large	with individual holders for the endoscope suitable for PCNL optic	
PD-FS-4001	ocular fix-focus black	modular ocular fix focus black with pin for easy mounting at the 3-joint-arm best for use with: 10.000 and 30.000 pixel optic	
PD-DS-1210	optic shifter	for all handles length 26mm	
PD-TD-0020	cleaning brush shifter set of 5	cleaning brush with eye for optic shifter PD-DS-1210 set of 5 pieces	
PD-PN-0015	All-Seeing Needle	consisting of: Needle, Working Shaft 4.8F, Mandrin disposable / sterile packed	
PD-PN-0005	8F Working Shaft	OD 2,66mm 8F ID 2,3mm 7F including Dilator autoclavable	
PD-PN-0020	10,5F Working Shaft	OD 3,46mm 10,5F ID 3,00mm 8,5F including Dilator autoclavable	
PD-PN-0010	Lunderquist guide wire	Size 0.035 x 80cm	

Micro Perc PCNL

Order No.	Article Description	Additional Description	Picture
PD-PN-0007	complete adapter set PCNL	consisting of: a) T-piece b) 3-port with female Luer-Lock c) Luer-Lock cap d) Tuohy Borst e) perfusion tube	
PD-IP-0100	Irrigation Pump IP 200	100ml/min 220 Volt incl. 5 pcs. PD-IP-0120 and 5 pcs. PD-IP-0121	
PD-IP-0120	Tubes for Irrigation Pump	1 piece	
PD-IP-0121	Spike for Tube	1 piece	